

BLACK PANTHER



REF	PM2810	PM2820	P2826
Ø / L (in mm)	4 / 13	11 / 2	26 / 2,5
Schaft / shank	HP	HP	HP
	5 / 1	5 / 1	1

BLACK PANTHER



REF	PM2910	PM2920	P2926
Ø / L (in mm)	4 / 13	11 / 2	26 / 2,5
Schaft / shank	HP	HP	HP
	5 / 1	5 / 1	1

Zirkon / Zircon

ZIRKO-SHINE

Diamantkörnung: grob
Diamond grit: coarse

1. Stufe
Abtragen und Formen

1. Step
Removing
and shaping

5.000 - 10.000 min⁻¹



REF	PM2510	PM2520	P2521
Ø / L (in mm)	4 / 13	11 / 2	17 / 2,5
Schaft / shank	HP	HP	HP
	5 / 1	5 / 1	5 / 1

ZIRKO-SHINE

Diamantkörnung: medium
Diamond grit: medium

2. Stufe
Glätten

2. Step
Smoothing

5.000 - 10.000 min⁻¹



REF	PM2610	PM2620	P2621
Ø / L (in mm)	4 / 13	11 / 2	17 / 2,5
Schaft / shank	HP	HP	HP
	5 / 1	5 / 1	5 / 1

High-End Premium Instrumente

für monolithisches oder teilverblendetes Zirkon, Lithium-Disilikat (LS2), Glaskeramiken und konventionelle Verblendkeramiken

- Zwei-Stufen-Finish zum Glanz
- Perfekte Oberfläche für den Antagonistenschutz
- Lange Standzeit (kein Trägerkern)
- Abgestimmte Formenauswahl für alle Bereiche
- sicher und effizient
- splitterfreie Bearbeitung scharfer Ränder und feinsten Bereiche bei Teilkronen, Inlays und Veneers

High-End premium instruments

for monolithic or partially veneered zirconium, lithium disilicat, hybrid-ceramics and all commonly used veneering porcelains

- two-step-finish for perfect shine
- perfect surfaces
- longlife instruments
- fine tuned range of shapes for all areas
- safe and efficient
- shatterproof preparation of edging and finest areas at crowns, inlays and veneers

Anwendungshinweis

Spezielles dreistufiges Poliersystem für Zirkonoxid mit erhöhter Standzeit. Für kleine nachträgliche Formkorrekturen mit anschließender Hochglanzpolitur.

Recommendations for use

Special 3-step polishing system for zirconium oxide with increased service life. For little subsequent shape corrections with following high-shine polishing.

ZIRKO-SHINE

Diamantkörnung: fein
Diamond grit: fine

3. Stufe
Hochglanzpolieren

3. Step
High-shine polishing

10.000 - 20.000 min⁻¹



REF	PM2710	PM2720	P2721
Ø / L (in mm)	4 / 13	11 / 2	17 / 2,5
Schaft / shank	HP	HP	HP
	5 / 1	5 / 1	5 / 1