

368 (257)

REF	368	368
Ø (1/10 mm)	018	023
Körnung / Grit		
	6	6
	FGXL	XL
L (mm)	5	5

379 (277)

REF	379	379	379
Ø (1/10 mm)	023	018	023
Körnung / Grit			
	6	6	6
	FGL	L	
	FGXL		XL
L (mm)	4,5	4,5	4,5

801 (001)

REF	801	801	801
Ø (1/10 mm)	010	012	014
Körnung / Grit			
	6	6	6
	FGXL	XL	XL
L (mm)	1,0	1,2	1,4

830L (234)

REF	830L
Ø (1/10 mm)	014
Körnung / Grit	
	6
	FGL
	L
L (mm)	5,0

879 (290)

REF	879
Ø (1/10 mm)	014
Körnung / Grit	
	6
	FGXL
	XL
L (mm)	10

Schaftarten / Shank Types

	FGL	FG-Schaft, lang / FG shank, long
	FGXL	FG-Schaft, extra lang / FG shank, extra long

Körnungen / Types of grit

	SC	Supergrob (180 µm) Supercoarse
	C	Grob (151 µm) Coarse
	M	Mittel (107-126 µm) Medium
	F	Fein (40 µm) Fine
	SF	Superfein (20 µm) Superfine
	UF	Ultrafein (10 µm) Ultrafine